


## Eliminate sample heating during common NMR measurements

Application Note 7



Nuclear Magnetic Resonance (NMR) measurements such as  $T_2$  distributions,  $T_1$ - $T_2$  and  $T_2$ -diffusion maps are employed to derive porosity, wettability, fluid types and other properties. These measurements are derived using the Carr-Purcell-Meiboom-Gill (CPMG) NMR pulse sequence [1]. This pulse sequence employs a train of radio frequency pulses which probe the exponential decay of the NMR signal. In order to probe the fastest decays, the pulses at the beginning of the CPMG pulse train must be close together. For slower decays, the pulse train must be long. The large number of radio frequency pulses will heat the sample through inductive heating which will alter the NMR decay itself. By smartly spacing the RF pulses, all the decays can be measured with far fewer pulses thus significantly reducing the heating.

In a typical CPMG sequence, the pulses in the train are linearly spaced. This results in trains with many closely spaced pulses, leading to heating of the sample. Heating is an issue because as the sample's temperature increases, its NMR signal decreases. This can lead to inaccuracies in the porosity derived from the NMR data. In addition, the petrophysical property being measured may be temperature dependent and thus there is no control of an important property. For porosity, we have explored the degree of inaccuracy and have found that heating can lead to changes in retrieved NMR porosity of approximately 10% in both bulk fluid and core samples. We have also explored the effect of salinity of the sample fluid.

One method to eliminate the unwanted heating of samples is to reduce the number of RF pulses by logarithmically spacing the echoes (Figure 1). This logarithmically spaced train will maintain the closely spaced pulses at the beginning of the sequence necessary to probe the smallest pores, while


decreasing the density of the pulses as the train gets longer. This leads to a reduction in the total number of pulses and therefore a reduction in the heating of the sample. We have implemented these logarithmically spaced CPMG sequences and will show that these sequences still accurately reproduce the NMR measurements for all samples, while eliminating sample heating.

### Key Findings:

- Traditional CPMG sequences can heat samples leading to errors in the porosity of up to 10%.
- The NMR response often includes a temperature dependency which will lead to difficulties in interpretation if the temperature is changed by the measurement itself.
- Log spaced CPMG eliminates sample heating by reducing number of pulses by up to 2 orders of magnitude.
- Log spaced CPMG measures the same decay response with far greater accuracy (porosity variation < 1%) and no temperature affects.



# Eliminate sample heating during high duty cycle NMR measurements



**Figure 1:** Linear vs. log spaced CPMG sequences. In the linear sequence the spacing between each 180 degree RF pulse is equal ( $\tau$ ). In the logarithmic spaced sequence the spacing between each RF pulse grows exponentially.

## Experiments:

All the NMR experiments were carried out using an Oxford Instruments GeoSpec 2+75 rock core analyser [2]. The pulse sequence development and data acquisition were performed using Green Imaging Technologies App Builder software [3]. Three sets of experiments were carried out to explore reducing heating during NMR measurements by employing a log spaced CPMG sequence. The first set of experiments explored the ability of the log spaced pulse sequence to accurately reproduce pore size distributions. The second set of experiments were designed to calibrate the change in NMR volume observed due to sample heating by the CPMG to the actual change in temperature of the sample. The final set of experiments explored the effect of heating of different types of samples (bulk vs. core and 2% vs. 20% brine saturation) via CPMG pulse sequences and examined to what degree log spaced sequences can reduce this heating.


To properly quantify the effects of heating by the CPMG sequences in experiments two and three, the volumes retrieved from the  $T_2$  distributions needed to be plotted as a function of time. To do this, new pulse sequences were developed where the CPMG  $T_2$  distribution measurement was put inside a loop. Each time through the loop, the CPMG pulse sequence was run several times. This raw NMR decay data was then averaged and a  $T_2$  distribution created from this averaged data. The measured volume was then retrieved by summing the area under the  $T_2$  distribution. The volume was then saved to a file along with a date/time stamp for each volume measurement. This data was then used to plot volume versus time. Changes in retrieved volume over time were examined as the sample was heated by the CPMG sequence. The same looping pulse sequence was used to retrieve the data of interest with both the log and linear CPMG pulse sequences.



### Experiment 1: Pore size distribution correlation

To test the ability of the log spaced CPMG sequence to accurately reproduce pore size distributions, three glass filled samples (2.5 cm diameter x 1 cm high) were stacked into the magnet. Each sample contained approximately 2.5 ml of 2% KCl brine. Each sample was doped with different amounts of copper sulfate which is paramagnetic and causes the principal  $T_2$  component for each sample's distribution to be shifted. The higher the copper sulfate concentration, the lower the peak  $T_2$  value. The three samples employed had peak  $T_2$  values at approximately 0.5 ms, 10 ms and 100 ms. When the three samples were stacked together in the magnet, they produce one combined  $T_2$  distribution with three distinct peaks. This distribution was then probed with both the linear and log CPMG sequences.

Figure 2 shows the comparison of the  $T_2$  distribution retrieved for the linear CPMG sequence (Figure 2 – blue trace) to that retrieved from the log spaced CPMG (Figure 2 – red trace). Table 1 shows the NMR parameters employed for both the linear and log spaced CPMG sequences employed. As mentioned in the experimental section, this  $T_2$ -distribution was retrieved from three samples with three different peak  $T_2$  values stacked on top of each other in the magnet. Figure 2 shows that there is excellent agreement between the distributions derived from the linear CPMG sequence as compared to the log CPMG sequence. The total volume of the three samples derived from the log and linear spaced sequences differs by less than 1%.



**Figure 2:** Comparison of  $T_2$  distributions retrieved from the linear spaced CPMG (blue trace) and log spaced CPMG (red trace) pulse sequences retrieved from three samples with three different peak  $T_2$  values stacked on top of each other in the magnet.

## Eliminate sample heating during high duty cycle NMR measurements



	Log Spaced CPMG	Linear Spaced CPMG
Tau ( $\mu\text{s}$ )	Log spaced between 100 $\rightarrow$ 1380	100
Number of Echoes	512	2500
Max $T_2$ (ms)	100	100
Recycle Delay (s)	0.5	0.5
Number of Scans	16	16
P90 ( $\mu\text{s}$ )	11.22	11.22
P180 ( $\mu\text{s}$ )	22.42	22.42
NMR Resonance (MHz)	2.457	2.457

**Table 1:** NMR parameters for comparison of performance of log vs linear spaced CPMG




### Experiment 2: Calibration of NMR volume to temperature variations

The second set of experiments set out to quantify the effect of CPMG heating by tracking the volumes measured as a function of time. As mentioned earlier, as a sample heats there will be a decrease in the volume measured from the NMR data. The blue dots in the upper panel of Figure 3 show the associated change in NMR volume as a bulk 2% KCl sample is probed as a function of time. The volume clearly drops from approximately 17.3 ml to 16.7 ml during the experiment. It should be noted that the scatter in the blue dots is associated with instability in the NMR data and is not truly caused by variations in temperature. The NMR volume data can also be plotted as a percentage change (Figure 3 – lower panel – blue dots). Also shown in the upper panel of Figure 3 is the sample temperature as a function of time. This data was recorded simultaneously with the NMR volume using an infrared thermometer and shows that the sample increases in temperature from approximately 27.50 C to 34.75 C. As expected, there was clearly an inverse correlation between the temperature of the sample and the observed NMR signal or volume. Specifically, the NMR signal is tied to temperature via a Boltzmann distribution. Therefore, the observed sample temperatures can be employed to predict the change in observed NMR volume.



This was completed for the temperature data recorded and the predicted percent change in NMR volume is plotted as the red dots in the lower panel of Figure 3. As can be seen there is good agreement between the measured change in NMR volume and predicted change in NMR volume. Any difference between the predicted and measured change in NMR volume can be attributed to the fact that the IR thermometer measures the surface temperature of the brine whereas the NMR measurement probes the entire bulk volume. The surface temperature could be cooler than the average temperature of the bulk sample as it is cooled by the air above the sample.




**Figure 3:** NMR volume of a bulk 2% KCl sample (blue dots – upper panel) measured as a function of time by a linear CPMG sequence. The temperature of the bulk sample (red dots – upper panel) was also logged simultaneously to the NMR measurement with an infrared thermometer. The lower panel shows the NMR volume plotted as a percent change (blue dots). The red dots in the lower panel are the percent change in NMR volume as predicted via the Boltzmann distribution using the temperature of the bulk sample.

# Eliminate sample heating during high duty cycle NMR measurements



## Experiment 3: Log versus linear CPMG heating changes

Once we were confident that the observed changes in NMR volume accurately reflected heating of the samples by the CPMG sequence, we began to study sample heating with both the log and linear CPMG sequences on various sample types. Table 2 summarizes the NMR parameters employed for the log and linear CPMG pulse sequences employed.



**Figure 4:** Plots of percent change in NMR volume measured with both the linear spaced CPMG pulse sequence (blue trace) and the log spaced pulse sequence (red trace). The log spaced CPMG sequence clearly reduces heating for both bulk samples (upper two panels) and saturated sandstone samples (lower two panels).

Figure 4 compares the percent change in NMR volume measured with the log and linear CPMG pulse sequences for the samples tested (Bulk 2% KCl, Bulk 20% KCl, sandstone saturated with 2% KCl, sandstone saturated with 20% KCl). For each of the samples, there is clearly more change in the volume measured for the linear CPMG sequence as compared to the log sequence. For example, for the upper right panel of Figure 4 the percent change in the measured NMR volume is approximately 7% for the linear CPMG sequence (blue dots). Conversely a percent change of only approximately 1% for the log CPMG sequence (red dots) is observed.

It should be noted that there remains a small change in the NMR volume measured by the log CPMG sequence. This indicates that despite having fewer echoes than the linear CPMG sequence, the log sequence still heats the sample to a certain degree. This heating however amounts to only approximately 1% change in the NMR volume measured.

The linearly spaced CPMG data presented in Figure 4 also allows interesting conclusions to be drawn on the difference in sample heating between bulk vs core samples and the effect of salinity on sample heating. Firstly, if we compare the linearly spaced CPMG data (blue dots) in the upper two panels of the



figure, we see that the higher salinity bulk sample heats more than the lower salinity sample (7% reduction in NMR volume for higher salinity sample vs. 4% reduction in NMR volume for lower salinity sample). This is expected because the increased salt concentration of 20% KCl brine lowers the specific heat capacity of the water allowing it to reach a higher temperature as compared to the 2% KCl brine. The same comparison of salinity for the core samples (Figure 4 – lower panels – blue dots) shows that again when saturated with the higher salinity brine the core sample heats more than when it saturated with the lower salinity brine. (8% reduction in NMR volume for higher salinity sample vs. 7% reduction in NMR volume for lower salinity sample). The effect however is not as pronounced as with the bulk samples.

	Bulk KCl Samples (2% and 20%)		KCl Saturated (2% and 20%) Sandstone Sample	
CPMG	Log	Linear	Log	Linear
Tau (μs)	Log spaced between 50 → 1795	50	Log spaced between 50 → 1795	50
# of Echoes	512	5000	512	5000
Max T <sub>2</sub> (ms)	100	100	100	100
Recycle Delay(s)	0.75	0.75	0.75	0.75
# of Scans	16	16	160	160
# of Meas.	1000	1000	100	100

**Table 2:** NMR parameters for study of heating by CPMG sequences in pore size distribution measurements.



# Eliminate sample heating during high duty cycle NMR measurements

Application Note 7



## Conclusion:

This paper presents results showing that the linearly spaced RF pulses in CPMG sequences (employed in  $T_2$  distributions and  $T_1$ - $T_2$  maps) can deposit excess energy into samples leading to unwanted heating. This heating can lead to inaccuracies in the measurements retrieved from the NMR data. This work has shown that depending on the type of sample (bulk vs sandstone core), the pulse parameters employed and the salinity of the brine, heating can lead to a discrepancy of up to 10% in the pore volume or porosity retrieved from  $T_2$  distributions.

A new pulse sequence was presented which eliminates the heating of samples by the many pulses in a linearly spaced CPMG pulse sequence. This new pulse sequence employs logarithmically spaced pulses. This new sequence maintains the closely spaced pulses at the beginning of the sequence necessary to probe the smallest pores while decreasing the density of the pulses as the train gets longer. This leads to a reduction (by as much as two orders of magnitude) in the total number of pulses necessary to accurately probe the exponential decay of the NMR signal and therefore a reduction in the heating of the sample. This work has shown that this new sequence can accurately reproduce  $T_2$  distributions and retrieve correct pore volumes from these distributions. The work also shows that the log spaced pulses in this modified CPMG sequence lead to no unwanted heating of the sample.

## References

1. Meiboom, S. and Gill, D., "Modified Spin Echo Method for Measuring Nuclear Relaxation Times", *Review of Scientific Instruments* (1958), **29**, 688-691
2. *Geo-Spec 2-75 User Manual, Version 1.8, Oxford Instruments.*
3. *GIT App Builder User Manual, Version 1.1 Green Imaging Technologies.*

## For further information:

M. Dick, D. Green, D. Veselinovic, "Reducing Sample Heating During NMR Measurements", *Society of Core Analysts Annual Symposium, Pau France, 2019*

M. Dick, D. Green, D. Veselinovic, "Accurate Pore Size Measurement Via NMR On Unconventionals", *Unconventional Resources Technology Conference, Colorado USA, 2019*

## Oxford Instruments Magnetic Resonance

For more information: [magres@oxinst.com](mailto:magres@oxinst.com)  
[www.oxinst.com/geospec](http://www.oxinst.com/geospec)

### UK

Tubney Woods, Abingdon, Oxfordshire, OX13 5QX, UK  
**Tel:** +44 (0) 1865 393 200 **Fax:** +44 (0) 1865 393 333

### USA

300 Baker Avenue, Suite 150, Concord, MA, 01742, USA  
**Tel:** +1 978 369 9933 **Fax:** +1 978 369 8287

### China

Floor 1, Building 60, No.461, Hongcao Road,  
Shanghai, 200233, China

**Tel:** +86 21 6073 2925 **Fax:** +86 21 6360 8535

## Green Imaging Technologies

For more information: [info@greenimaging.com](mailto:info@greenimaging.com)  
[www.greenimaging.com](http://www.greenimaging.com)

### Canada

520 Brookside Drive, Suite B, Fredericton, NB, E3A 8V2,  
Canada

**Toll Free:** +1 888 944 8462

**Tel:** +1 506 458 9992 **Fax:** +1 506 458 9615

This publication is the copyright of Oxford Instruments plc and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trademarks and registrations. © Oxford Instruments plc, 2019. All rights reserved. Ref: MR/211/0619



The Business of Science®